



# Mobile Independence Pty Ltd

ABN: 90 531 747 431



3/16 Commercial Drive  
Caboolture, QLD 4510

Ph: 07 3881 3868

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: [reception@mobileindependence.com.au](mailto:reception@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Air Comfort Chair



The ultimate in comfort and mobility, the Comfort Chair is an easy care solution for elderly or special needs seating. Cushioned adjustable air cell support has been incorporated into the design. The four lateral support wings fold away for easy transferring of the patient. This chair has a reclinable seat along with an adjustable backrest and calf rest, tray, neck support and footplate. The frame is completely constructed from sturdy stainless steel.

### Features:

- Adjustable backrest (90° - 170°) and calf rest
- Reclinable seat
- Air system for pressure relief
- Full stainless steel frame
- Three lockable castors and one steering castor (non-marking)
- Breathable polyurethane cover (fire retardant)
- Replaceable and easy to clean cover
- Four folding wings
- Standard dining tray
- Magazine pocket

### Specifications:

Seat Height:	63cm
Seat Width:	51cm
Seat Depth:	50cm
Overall Width:	73cm
Overall Length:	135 - 188cm
Table to Floor:	88cm
Table Slide:	23cm
Table Size:	40 x 40.5 x 2cm
Back Reclining:	102° - 136°
Seat Tilting:	0 - 18°
Leg Elevating:	70° down to 5° up
Weight of Chair:	43.5kg
Maximum User Weight:	120kg

10 years warranty on frame