

Air Comfort Chair



The ultimate in comfort and mobility, the Comfort Chair is an easy care solution for elderly or special needs seating. Cushioned adjustable air cell support has been incorporated into the design. The four lateral support wings fold away for easy transferring of the patient. This chair has a reclinable seat along with an adjustable backrest and calf rest, tray, neck support and footplate. The frame is completely constructed from sturdy stainless steel.

Features:

Adjustable backrest (900-1700) and calf rest

Reclinable seat

Air system for pressure relief

Full stainless steel frame

Three lockable castors and one steering castor (non-marking)

Breathable polyurethane cover (fire retardant)

Replaceable and easy to clean cover

Four folding wings


Standard dining tray

Magazine pocket

10 years warranty on frame



Unit 1 / 60 Kremzow Road Brendale Qld 4500 www.mobileindependence.com.au

 (07) 3881 3868 Fax (07) 3881 0525 reception@mobileindependence.com.au



Mobile Independence