



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

All Day Slings with Head Support



A very comfortable and fully supportive sling designed for standard transfers including from the floor. The mesh material is suitable for bathing.

Specifications

	Nylon Code	Mesh Code	A	B	C
Small	HC01-1500S	HC01-1550S	880mm	1350mm	750mm
Medium	HC01-1500M	HC01-1550M	950mm	1400mm	800mm
Large	HC01-1500L	HC01-1550L	1050mm	1550mm	900mm
Extra Large	HC01-1500XL	HC01-1550XL	750mm	800mm	900mm

Warranty: 1 year limited