



# Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## All Day Slings



Similar to the cradle sling without head support and backing foam. The patient requires some upper body control to support their head.

### Specifications

	Nylon Code	Mesh Code	A	B	C
<b>Small</b>	HC01-1300S	HC01-1350S	500mm	1000mm	750mm
<b>Medium</b>	HC01-1300M	HC01-1350M	550mm	1030mm	800mm
<b>Large</b>	HC01-1300L	HC01-1350L	610mm	1150mm	900mm
<b>Extra Large</b>	HC01-1300XL	HC01-1350XL	650mm	1180mm	1000mm

**Warranty:** 1 year limited