

# Jazzy<sup>®</sup>

## SELECT SERIES

Jazzy Select<sup>®</sup> Power Chairs offer the advantages of inline motor technology for enhanced efficiency, torque and range. The highly manoeuvrable Jazzy Select is loaded with performance and convenience features, while the Jazzy Select 6 provides maximum front stability for superior performance indoors and out. The Jazzy Select 14 XL completes the package, with a powerful 60 amp PG VR2 controller.

**Jazzy Select<sup>®</sup> 6**  
shown in Jazzy Red



**Jazzy Select<sup>®</sup>**  
shown in Jazzy Red



**Jazzy Select<sup>®</sup> 14 XL**  
shown in Jazzy Red



21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • [www.pridemobility.com](http://www.pridemobility.com)

**Pride<sup>®</sup>**  
Mobility Products  
Australia Pty. Ltd.

SPECIFICATIONS	Jazzy Select®	Jazzy Select® 14XL	Jazzy Select® 6
<b>Weight capacity</b>	136 kg (300 lbs.)	181 kg (400 lbs.)	136 kg (300 lbs.)
<b>Maximum speed</b>	Up to 6 km/h (4 mph)	Up to 6.5 km/h (4mph)	Up to 6 km/h (4 mph)
<b>Ground clearance</b>	4.76 cm (1.875")	7.6 cm (3")	5 cm (2")
<b>Turning radius</b>	51 cm (20")	53 cm (21")	52 cm (20.5")
<b>Overall length<sup>1</sup></b>	90 cm (35.5")	96.5 cm (38")	93 cm (36.625")
<b>Overall width</b>	60 cm (23.75")	63 cm (24.75")	60 cm (23.5")
<b>Wheels</b>	Drive: 25 cm (10")	Drive: 35.5 cm (14")	Drive: 25 cm (10") solid
<b>Front castors</b>	13 cm (5") solid anti-tips	15 cm (6") solid anti-tips	15 cm (6") solid
<b>Rear castors</b>	15 cm (6") solid	20 cm (8") solid	15 cm (6") solid
<b>Suspension</b>	Active-Trac®	Limited	Active-Trac®
<b>Drive train</b>	Two-motor, in-line, mid-wheel drive	Two-motor, mid-wheel drive	Two-motor, in-line, mid-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)		
<b>Standard electronics</b>	40 amp, PG GC 2	60 amp, PG VR2	40 amp, PG GC 2
<b>Specialty controls</b>	N/A		
<b>Battery charger</b>	Off-board, 3 amp	Off-board, 4-amp	Off-board, 3 amp
<b>Per-charge range</b>	Up to 24 km (15 miles)		
<b>Battery requirements</b>	(2) 12 volt, deep cycle		
<b>Battery size</b>	35 AH	35 AH	35 AH
<b>Battery weight</b>	11 kg (24.5 lbs.) each	11 kg (24.5 lbs.) each	11 kg (24.5 lbs.) each
<b>Base weight</b>	45 kg (99 lbs.)	48 kg (106.5 lbs.)	45 kg (100 lbs.)
<b>Standard seat weight</b>	19 kg (42.5 lbs.) (high-back)	19 kg (42.5 lbs.) (high-back)	19 kg (42.5 lbs.) (high-back)
<b>Max. seat dimensions</b>	46 cm x 46 cm (18" x 18")	46 cm x 51 cm (18" x 20")	46 cm x 51 cm (18" x 20")
<b>Warranty</b>	2-year limited		

## COLOUR CHOICE



Jazzy Red

1) Without front riggings



## FEATURES (Jazzy Select®)

- Dual in-line motors for enhanced efficiency, torque, range and performance
- New 40 amp, PG GC 2 controller
- Build with ease of service in mind
- Easy side access to batteries
- Tool-less front panel access to electronics to minimize service time

## FEATURES (Jazzy Select® 14 XL)

- Patented mid-wheel drive design for tight space manoeuvrability
- 35.5 cm (14") knobby tyres and optimum weight distribution provide superior traction
- Upgradeable to 55 AH batteries for maximum range per charge
- New 60 amp, PG VR2 controller
- Off-board, 4 amp charger

## FEATURES (Jazzy Select® 6)

- 15 cm (6") casters on front and rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear caster beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform
- Built with ease of service in mind

## OPTIONS

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder

21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; due to the Pride's commitment to constant improvement and development, We reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

**PRIDE**<sup>®</sup>  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST™ —