



# Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868  
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au) Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Aluminium Shower Chair



This shower is a light weight, height adjustable, aluminium frame designed to assist patients in showering themselves.

### Features

Height adjustable  
Back rest optional

### Specifications

<b>Weight Capacity</b>	105 kg
<b>Height:</b>	47 cm – 59.5 cm
<b>Width:</b>	45.5 cm
<b>Length:</b>	33 cm
<b>Weight:</b>	2.8 Kg

**Warranty:** 1 year limited