



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Bariatric Walker



This model is designed for the larger person, with a wider frame and greater weight capacity than a standard Walker. It has the cable brake system which operates by applying a light squeeze on the brake lever to temporarily stop, or by pushing down on the brake lever to lock the brakes before sitting on the comfortable padded seat. Supplied with a detachable backrest, the frame folds up to make it convenient for transport and storage. The steel frame weighs 12.10 kg and can support someone up to 181kg in weight.

Features

- Fold up frame for ease of transport and storage
- Steel construction with powder coated finish
- Detachable height-adjustable handles
- Foam covered detachable backrest
- Wire shopping basket sits under seat

Specifications

Weight Capacity	181 kg
Seat height:	55 cm
Width:	77 cm
Height:	81.7-97.2 cm
Length:	70 cm
Seat Width:	46 cm
Seat Depth:	35 cm
Weight:	12.1 kg

Colour: Blue only

Warranty: 1 year limited