

Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525


QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Bed Pole



The Bed Pole is designed to sit under the mattress and provide a firm handle to assist with moving and positioning in bed.

Features

Made of a powder coated white 1" steel tubing
Easy to install, simply slide the "s" shape base between the mattress and the bed base

Specifications

Height:	50.5 cm
Width:	42 cm
Length:	38 cm
Weight:	1.8 Kg

Warranty: 1 year limited