



# Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868  
Fax: 07 3881 0525




QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Bed Rail



This Bed Rail is designed to sit under the mattress to provide a firm support for moving and positioning within the bed. It can also provide a level of security and comfort for some people and reduce the risk of falling out of bed.

### Features

Made of a powder coated white 1" steel tubing  
Suitable for most beds  
Easy to install, simply slide between the mattress and the bed base

### Specifications

**Height:** 48.5 cm  
**Width:** 60 cm  
**Length:** 58 cm  
**Weight:** 2 Kg

**Warranty:** 1 year limited