



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Boomer Lifter



The Boomer is our highest capacity bariatric hoist with a Safe Working Load of 350Kg. It has an extra wide base for stability and wide handles to accommodate 2 operators. It is fitted with a four point hanger designed for use with the Boomer bariatric slings.

Features

It is fitted with a four point hanger designed for use with the Boomer bariatric slings
Haycomp hoists are designed to pick up from the floor and lift to clear 1 metre.

Specifications

Weight Capacity	350 kg
Lift:	Electric
Leg opening:	Electric
Front leg height:	15.5 cm
Base width:	98 cm
Overall length:	141 cm
Mast height:	140 cm
Floor lift:	Yes

Options

Weigh scale

Warranty: 2 year limited