



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Cradle Sling (Toileting)



The sling is cut away at the base to allow access for toileting. Extra support is provided by soft under arm pads and a chest strap.

Options

Also available without head support

Specifications

	Nylon Code	Mesh Code	A	B	C
Small	HC01-1200S	HC01-1250S	720mm	1250mm	900mm
Medium	HC01-1200M	HC01-1250M	750mm	1330mm	1000mm
Large	HC01-1200L	HC01-1250L	800mm	1450mm	1120mm
Extra Large	HC01-1200XL	HC01-1250XL	850mm	1550mm	1200mm

Warranty: 1 year limited