



# Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868

Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: [reception@mobileindependence.com.au](mailto:reception@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Easy Grip Utensils



### Features

- The spoon, fork and knife have a large easy to grip handle which can be adapted for left or right hand use
- The spoon, fork and knife head can be removed and rotated to allow the user to position the spoon in the most comfortable hand position
- Made of stainless steel with a moulded plastic hand grip
- Dishwasher safe

### Specifications

<b>Overall Length:</b>	210 mm
<b>Handle Length:</b>	120 mm
<b>Knife Weight :</b>	75 g
<b>Fork Weight :</b>	70 g
<b>Spoon Weight:</b>	80 g

**Colour:** Black

**Warranty:** 1 year limited