



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Economy Walker 8" Wheels



4 wheeled, Aluminium, Loop Brake, 8" Wheel.

Features

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 8" x 1" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard

Specifications

Weight Capacity	136 kg
Seat height:	58.4 cm
Width:	62.5 cm
Height:	81.7-97.2 cm
Length:	65.4 cm
Seat Width:	30.5 cm
Seat Depth:	3.5 cm
Weight:	6 kg

Colour: Red

Warranty: 1 year limited