



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Emu B



The Emu B is our highest capacity bariatric stand hoist with a Safe Working Load of 400Kg. It has an extra wide base for stability and wide handles to accommodate 2 operators..

Features

Specifications

Weight Capacity	400 kg
Lift:	Electric
Leg opening:	Electric
Front leg height:	15.5 cm
Base width:	96 cm
Overall length:	139 cm
Mast height:	151 cm
Floor lift:	No

Warranty: 2 year limited