



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Emu Walker



HAYCOMP produce a product that is designed to reduce the injury to nurses and persons required to move patients with ease..

Features

- It has an extra deep knee pad to provide maximum comfort and the elimination of foot slippage.
- It has high lift and a unique handle system for greater comfort and safety.
- The Emu features an ergonomic design, unique slings and accessories to assist in high quality patient care.
- Australian Made

Specifications

Weight Capacity	200 Kg
Lift:	Electric
Leg Opening:	Electric
Front Leg Height:	14 cm
Base Width:	65 cm
Overall Length:	110 cm
Mast Height:	147 cm
Floor Lift:	Yes

Accessories

Knee strap
Slings
Bum Strap
Rehab Arms

Warranty: 2 year limited