



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Folding Canes



These are height adjustable Folding Canes fold into 4 convenient parts making it easy to take with you wherever you go. The cane folds to 29cm in length when collapsed.

Features

- Supplied with carrying case
- Automatically opens to a full size cane
- Height adjustable in 25mm (1") increments
- Constructed with 22mm anodised aluminium
- Secured with locking ring for added safety
- Non slip rubber tips (with metal inserts) ensure a secure grip

Specifications

Weight Capacity	114 kg
Height:	76 - 99 cm
Weight:	370 g

Colour: Black, Bronze or Chrome

Warranty: 1 year limited