

# Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868  
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Forearm Walker Economy



The Gutter Arm Walker is a powder coated mild steel frame folding walker designed to support the user in an upright position to assist with mobility and rehabilitation. It has four independent swivelling castors and can be easily folded for transport and storage. The arm supports are height adjustable and the handgrips can be moved forward or backward to suit the individual.

### Features

- Lightweight steel frame
- Easy to fold, transport and store
- Padded vinyl covered forearm supports
- Four independent swivelling castors
- Height adjustable
- Handgrips are adjustable

### Specifications

<b>Weight Capacity</b>	114 kg
<b>Height:</b>	95-130 cm
<b>Length</b>	80 cm
<b>Width</b>	63 cm
<b>Weight</b>	12.2 kg

**Colour:** Grey

**Warranty:** 1 year limited