

Mobile Independence Pty Ltd

ABN: 90 531 747 431



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Forearm Walker



4 wheeled, Steel, Loop Brake, 8" Wheel, Red/Black Patterned.

Features

- Steel frame
- Easy to Fold
- Four 8" x 2" wheels
- Padded Gutter arms
- Loop brake is ideal for slow down and permanent stop
- Wire Basket
- Clear acrylic tray
- Adjustable brake lever height

Specifications

Weight Capacity	120 kg
Seat height:	60 cm
Width:	58 cm
Height:	121 cm
Length:	68 cm
Seat Width:	37.5 cm
Seat Depth:	37.5 cm
Weight:	14 kg

Colour: Red/Black Patterned

Warranty: 1 year limited