

GOCHAIR TRAVEL CHAIR

The first travel chair that combines super portability, manoeuvrability and style at a low, low cost.



Specifications	
Drive wheels	20 cm (8") solid
Front casters	13 cm (5") solid
Rear anti-tips	8 cm (3") (solid)
Maximum speed	5.5 km/h (3.5 mph)
Ground clearance	4 cm (1.75")
Turning radius	57 cm (22.5")
Length	76 cm (30")
Width	47 cm (18.5")
Seat type	Foldable compact lightweight seat
Battery requirements/size	(2) 12V, 12 amp
Range	Up to 16 km (10 miles)
Charger	Off-board, 2 amp
Electronics	50 amp, PG VSI
Weight capacity	113 kg (250 lbs.)
Base weight	27 kg (59.5 lbs.) (front and rear assembled without batteries)
Standard seat weight	13.5 kg (30 lbs.)
Battery pack weight	9 kg (20 lbs.)
Suspension	Limited
Heaviest Piece (Rear)	16.5 kg (36.5 lbs.)
Warranty	2 year limited

PORTABILITY

- Easy to remove battery pack and compact, lightweight frame
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle boots

Available colours:



Red Blue

PERFORMANCE

- Top speed of 5.5 km/h (3.5 mph) and enough power to get you where you want to go
- Twin motors for traction and precision handling

MANOEUVRABILITY

- Tight turning radius of 57 cm (22.5")



"The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice".

"Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement".