

## SPECIFICATIONS

	Go-Go® Ultra X	Go-Go® Elite	Go-Go® Elite Plus	Go-Chair®
Max. weight capacity:	260 lbs.	275 lbs.	300 lbs./325 lbs. w/HD version	250 lbs.
Maximum speed: <sup>1</sup>	Up to 4 mph	Up to 4.25 mph	Up to 4.5 mph	Up to 3.5 mph
Ground clearance: <sup>2</sup>	3W: 1.375" / 4W: 1.375"	2.5"	3"	1.75"
Turning radius: <sup>2</sup>	3W: 32.5" / 4W: 44.125"	3W: 32.5" / 4W: 43.25"	3W: 34.5" / 4W: 45.75"	26.75" without foot platform
Overall length: <sup>2</sup>	3W: 37" / 4W: 39.75"	3W: 36" / 4W: 39.75"	3W: 40" / 4W: 40.25"	29" without foot platform
Overall width: <sup>2</sup>	3W: 19.5" / 4W: 19.5"	19.25"	3W: 21.125"	19"
Battery charger:	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp
Per-charge range: <sup>1,5</sup>	Up to 8 miles (w/ 2 x 12 AH)	Up to 8 miles (w/ 2 x 12 AH)	Up to 12 miles	Up to 8 miles
Compatible batteries:	12 AH	12 AH: 18 AH (optional)	18 AH	(2) 12V, 12 AMP
Weight:				
Batteries (2 required):	12 AH (20 lbs. battery pack) <sup>4</sup>	3W: 12 AH (20 lbs. battery pack) <sup>4</sup> 4W: 18 AH (33 lbs. battery pack) <sup>4</sup>	(2) 12V deep-cycle	(2) 12V, 12 AMP
Weight w/o battery pack: <sup>2</sup>	3-wheel: 72 lbs. 4-wheel: 78.5 lbs.	3-wheel: 75 lbs. 4-wheel: 78.5 lbs.	3-wheel: 84 lbs. 4-wheel: 87 lbs.	
Heaviest piece: <sup>2</sup>	3-wheel: 28 lbs. (front section) 4-wheel: 33.5 lbs. (front section)	3-wheel: 28 lbs. (front section) 4-wheel: 32.5 lbs. (front section)	3-wheel: 32.5 lbs. (rear section) 4-wheel: 35 lbs. (front section)	34.5 lbs. (rear section)

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup>Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

<sup>3</sup>AGM or gel cell type recommended.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications



## Dealer Information

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —

GOGOTRI\_07\_10

**GO GO**  
TRAVEL MOBILITY



**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —

**GO GO**<sup>®</sup>  
ULTRA X

**GO GO**<sup>®</sup>  
ELITE TRAVELLER<sup>®</sup>

**GO CHAIR**<sup>®</sup>

&

ELITE TRAVELLER<sup>®</sup> PLUS/HD



### FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride<sup>®</sup>'s exclusive black, non-scuffing tires
- Auto-connecting front to rear cable



### FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride<sup>®</sup>'s exclusive black, non-scuffing tires
- Auto-connecting front to rear cable
- Includes three sets of easily changeable colored panels in three exciting colors - Red, Blue and Silver



### PORTABILITY

- Easy to remove battery pack
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle trunks

### PERFORMANCE

- Top speed of 3.5 mph and enough power to get you where you want to go
- Twin motors for traction and precision handling

### MANEUVERABILITY

- Compact design for a tight turning radius