



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525




QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

General Purpose No Head Support Yoke Sling



For use with the WOMBAT, JOEY & JUMBUCK Patient Lifting Machines • Yoke Bar Attachment required • S.W.L 200kgs

Options

Also available with head support

Specifications

	Nylon Code	Mesh Code	A	B	C
Small	HC01-1000S	HC01-1050S	530mm	990mm	740mm
Medium	HC01-1000M	HC01-1050M	570mm	1170mm	840mm
Large	HC01-1000L	HC01-1050L	650mm	1350mm	870mm
Extra Large	HC01-1000XL	HC01-1050XL	680mm	1400mm	980mm

Warranty: 1 year limited