



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Hi Lo Bed



Designed to look like an ordinary mattress and base, with all the features of a hospital bed

Features

- Move off and on to bed with ease • Zero Gravity
- Adjustable level to assist care providers
- Whisper quiet, low voltage DC lift system
- Independent head and foot lift

Specifications

Weight Capacity	250 kg
Height	42 – 85 cm

Options

- Water proof health cover
- Aluminium rail set
- Pocket spring pillow
- Latex pillow
- Pocket spring mattress
- Choice of fabric
- Long single or King single

Warranty: 10 year Structural and 5 year Electrical