

# Mobile Independence Pty Ltd

ABN: 90 531 747 431



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

Email: [reception@mobileindependence.com.au](mailto:reception@mobileindependence.com.au)

## High Back Chair



Adjustable High Back Chair.

### Features

- High back chair
- Adjustable height
- Contoured back rest
- Padded arm rests
- Vinyl upholstery
- Powder coated steel frame

### Specifications

- Overall Width: 640mm
- Seat Width: 490mm
- Seat Depth: 480mm
- Seat Height: 455 - 605mm
- Back Height: 730mm

**Colour:** Teal and Mocha

**Warranty:** 1 year limited