



Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Instrument Trolley



Every Oxford trolley drawer is now fitted with full extension, roller bearing runners. Access the whole drawer with silky smooth operation. Drawers are full stainless steel or powder coated mild steel construction for ultimate durability.

Features

- Optional rails
- 2 shelves with 48cm clearance or 3 shelves with 30cm clearance
- Wheels: standard 7.5cm diameter
- Strong but lightweight welded construction
- Standard colour blue weave
- Manufactured from 304 grade stainless steel with satin finish

Specifications

Width:	50 cm
Height:	90 cm
Length:	50 cm
Shelves:	48 cm

Warranty: 1 year limited