

Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Joey Lifter



The Joey is a proven performer in institutional lifting, with electric lift and leg opening. It's a very popular hoist for Hostel, Nursing Home and Hospital use.

Features

- Tilt cradle frame. This provides maximum comfort and safety for patient and carer.
- Standard 100mm castors are supplied as standard for easy rolling.
- Haycomp hoists are designed to pick up from the floor and lift to clear 1 metre. The clearance height may vary with different hanger and sling combinations.
- Australian Made

Specifications

Weight Capacity	200 Kg
Lift:	Electric
Leg Opening:	Electric
Front Leg Height:	14 cm
Base Width:	65 cm
Overall Length:	110 cm
Mast Height:	136 cm
Floor Lift:	Yes

Specifications

Power Tilt

It can be supplied with a yoke hanger for use with GP slings

Weigh scale

Smaller front castors can be supplied when under bed clearance is an issue.

Warranty: 2 year limited