



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Jumbuck Walking Sling



For use with the Jumbuck Patient Lifting Machine fitted with the Walking Accessories • S.W.L 240kgs

Specifications

	Nylon Code	A Waist Pad Height	B Waist Pad Length	Waist Size
Small	HC07-1800S	290mm	750mm	700-1050mm
Medium	HC07-1800M	290mm	920mm	800-1200mm
Large	HC07-1800L	290mm	1100mm	1000-1400mm
Extra Large	HC07-1800XL	290mm	1300mm	1200-1600mm

Warranty: 1 year limited