



Mobile Independence Pty Ltd

ABN: 90 531 747 431



3/16 Commercial Drive
Caboolture, QLD 4510

Ph: 07 3881 3868

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Jumbuck 240 Lifter



The Jumbuck 240 is a higher capacity multipurpose hoist, designed for special care Nursing homes and Hospitals.

Features

- It has a selection of practical accessories that are easy to change.
- It has electric lift and leg opening
- Standard 100mm castors are supplied as standard for easy rolling. Smaller front castors can be supplied when under bed clearance is an issue.
- Haycomp hoists are designed to pick up from the floor and lift to clear 1 metre. The clearance height may vary with different hanger and sling combinations.
- Australian Made

Specifications

Weight Capacity	240 Kg
Lift:	Electric
Leg Opening:	Electric
Front Leg Height:	14 cm
Base Width:	65 cm
Overall Length:	125 cm
Mast Height:	165 cm
Floor Lift:	Yes

Options

Weigh scale

Tilt cradle frame. For use with slings that provide the maximum comfort and safety

Yoke bar hanger for use with GP slings

The stretcher frame attachment enables the lifting of Jordon frame and some Scoop stretchers

Rehab walking attachments. These include a walking hanger bar, walking sling and walking handles

Warranty: 2 year limited