

# Lightweight wheelchair L406-16

## Features

Attractive modern paint finish

Self-propel wheel chair with attendant brakes & folding back & quick release wheels

Swing away armrest, swing in swing out leg rests, total weight of 13kg

Washable anti-bacteria padded upholstery which is removable

Adjustable seat depth



Colour – Flame Red

## Technical Specs - Australian Standards Certified (AS/NZS 3695.1:2011)

Overall Width 63cm

Overall Length 107cm

Overall Height 92 cm

Seat Width 40cm

Seat Depth 38-43cm

Seat height from ground 48cm

Weight Capacity 115kg

Total weight 13kg