

Pride® Lift Chairs

END-TO-END

A two-position is a TV recline.
 A three-position is a full recline.
 An infinite-position is a full recline with independent legrest function.




Petite to Small Medium to Large Medium to Large Medium to Large Medium to Large Medium to Large Small to Medium Medium to Large Medium to Large Small to Medium

Pride Lift Chairs

The Ultimate in Comfort & Style!

	C1	C5	C6	C15	D30	550	560	585	805	660
	3-Position Lift Chair with fully removable upholstery components 163 cm & below (5'4" & below)	3-Position Chaise Lounger with fully removable upholstery components 163 cm (5'4") & below	Infinite Position Chaise Lounger with fully removable upholstery components 163 cm (5'4") & below	3-Position Lift Chair 157.5 cm - 178 cm (5'2" - 5'10")	3-Position Chaise Lounger 163 cm (5'4") & below	3-Position Chaise Lounger (Heat & Massage option now available) 163 cm - 183 cm (5'4" - 6'0")	3-Position Chaise Lounger 163 cm (5'4") & below	3-Position Chaise Lounger 165 cm - 183 cm (5'5" - 6'0")	2-Position Chaise Lounger Wall Hugger that requires only 10 cm (4") of space between the wall and chair 163 cm - 183 cm (5'4" - 6'0")	Infinite-Position Chaise Lounger 163 cm (5'4") & below
Weight Capacity	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300lbs.)	147 kg (325lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	136 kg (300lbs.)	170 kg (375 lbs.)
Back Width	51 cm (20")	71 cm (28")	71 cm (28")	71 cm (28")	74 cm (29")	89 cm (35")	71cm (28")	81 cm (32")	76 cm (30")	71 cm (28")
Seat to Floor	41 cm - 46 cm (16" - 18")	46-51cm (18-20")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	48cm (19")	51 cm (20")	48 cm (19")	48 cm (19")
Seat Depth	51 cm (20")	46-51 cm (18-20")	53 cm (21")	48 cm (19")	56 cm (22")	51 cm (20")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")
Seat Width	46 cm (18")	56 cm (22")	56 cm (22")	53 cm (21")	56 cm (22")	53 cm (21")	48 cm (19")	51 cm (20")	51 cm (20")	48 cm (19")
Top of Back to Seat	66 cm (26")	66-71 cm (26-28")	66 cm-71cm (26-28")	69 cm (27")	66 cm (22")	69 cm (27")	63.5 cm (25")	66 cm (26")	69 cm (27")	63.5 cm (25")
Back Style	Blown Fibre Split-Back	Blown Fibre Sewn Pillow-back	Blown Fibre Sewn Pillow-back	Blown Fibre Button-back	Blown Fibre Split-back	Blown Fibre Split-back	Blown Fibre Pillow Back	Blown Fibre Biscuit-Back	Blown Fibre Sewn Split-Back	Blown Fibre Pillow Back
Fully Padded Chaise	Yes ^{1, 2}	Yes ^{1, 2}	Yes ^{1, 2}	No	Yes ¹	Yes	Yes ¹	Yes	Yes	Yes ¹
Integrated Battery Backup³	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Recline Positions	3	3	Infinite	3	3	3	3	3	2	Infinite
Standard Arm Covers	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Standard Head Cover	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Colour Options	Arctic Blue, Chocolate, Oatmeal	Arctic Blue, Chocolate, Oatmeal, Levante (vinyl): Royal Blue	Arctic Blue, Chocolate, Oatmeal, Levante (vinyl): Royal Blue	Blue, Vino, Pine, Maple, Rivergum, Espresso, Levante (vinyl): Burgundy, Royal Blue, Pine Mountain, Mushroom	Blue, Vino, Pine, Rivergum, Levante (vinyl): Burgundy, Royal Blue, Pine Mountain	Almond, Rivergum, Merlot	Blue, Vino, Pine	Blue, Vino, Levante (vinyl): Mushroom	Blue, Espresso	Rivergum

FABRIC SELECTIONS



LEVANTE VINYL



Note: Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1. Not all fabrics are available for every model. Please see Specifications for details.