



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Medication Trolley



Oxford medication trolleys provide a range of solutions to almost any situation. Trolleys may be configured to a range of unit sizes from 12 through to 50.

Features

- Trolley with rails 3 sides
- Drawers: 3 x 10cm and 1 x 25cm deep
- 3 top drawers have removable aluminium compartments each with 3 dividers

Specifications

Width: 50 cm
Height: 100 cm
Length: 60 cm
Draws: 3 x 10 cm 1 x 25 cm

Warranty: 1 year limited