



Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Mobile Chair Shower Commode Drop Arm



The Commode combines a manoeuvrable and stable construction with a durable, stainless steel manufacture. Designed and built to the rigorous Australian Standard. All models feature washable surfaces, soft padded seat, PVC back rest, 12.5mm diameter stainless steel and reinforced Nylon castors. It is flexible in design, with user comfort and safety been a primary consideration. Forward location of front wheels ensures stability even when the user leans forward. Stable even when standing on the forward located foot plate.

Features

- Soft padded seat (open or closed front)
- Drop arm
- Pan holder
- Sliding footrest with extra grip surface
- Wheels: 12.5cm standard with 2 locking
- Very stable
- Australian Made

Specifications

Weight Capacity	120Kg
Width:	65cm
Height:	110 cm
Length:	70 cm
Material:	Stainless Steel

Warranty: 1 year limited