

Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Over Toilet Aid



The Over Toilet Aid is height adjustable and designed to assist patients in lowering themselves onto the toilet, and rising from it.

Features

- Zinc (anti corrosion) pre-treatment
- Height Adjustable 1" increments
- Supplied standard with seat, lid, splashguard and pail
- Powder coated white frame
- Can be used as a commode

Specifications

Weight Capacity	130 kg
Seat height:	48 cm - 61cm
Width:	53 cm
Length:	57 cm (22.5")
Seat Width:	48cm
Weight:	5.2 Kg

Warranty: 1 year limited