

Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Over Toilet Aid Folding



The Folding over Toilet Aid is height adjustable and designed to assist patients in lowering themselves onto the toilet, and raising from it. When not in use it easily folds for convenient storage.

Features

- Height adjustable
- Zinc pre-treatment
- One-piece blow moulded plastic seat
- Supplied with seat, pail and splashguard
- Micro adjustment for uneven floors
- Supplied fully assembled

Specifications

Weight Capacity	114 kg
Seat height:	51cm to 63cm
Width:	59 cm
Length:	49 Cm
Seat Width:	45 cm
Weight:	5.2 Kg

Warranty: 1 year limited