

Pride® Walkers

PMW 462

4-wheeled, Aluminium, Loop Brake, 6" Wheel

TECHNICAL SPECIFICATIONS

Seat height:	55 cm (21.5")
Width:	63.5 cm (25")
Height:	85-95 cm (33.5"-37.5")
Length:	57 cm (22.5")
Seath Width:	35 cm (14")
Seat Depth:	31 cm (12")
Weight:	6.5 kg (14.5 lb)
Weight Capacity:	114 kg (250 lb)
Colour:	Red/Black Patterned
Warranty:	1 year limited

FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 6" x 1" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard



PMW 463

4-wheeled, Aluminium, Loop Brake, 8" Wheel

TECHNICAL SPECIFICATIONS

Seat height:	55 cm (21.5")
Width:	63.5 cm (25")
Height:	88-98 cm (33.5"-38.5")
Length:	60 cm (23.5")
Seath Width:	35 cm (14")
Seat Depth:	31 cm (12")
Weight:	7 kg (15.5 lb)
Weight Capacity:	114 kg (250 lb)
Colour:	Red/Black Patterned
Warranty:	1 year limited

FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 8" x 1.5" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard



Pride®
Mobility Products
Australia Pty. Ltd.

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

21 Healey Rd.
Dandenong, Victoria 3175
www.pridemobility.com

ACN 088 609 661 ABN 66 088 609 661

Tel: 03 9706 4611
Fax: 03 9706 4622
www.pridemobility.com



ISO 9001:2008
SAI Global