



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Portable Bariatric Scale



Portable scale for the bariatric with two side entries.

Features

- 4 x 375kg load cells beneath the pan (pictured)
- The indicator is a Rice Lake 120PLUS with an internal rechargeable battery
- The 2 entry sides have a bevelled edge to prevent a tripping point
- Weigh plate powder coated
- The rails are stainless steel & very sturdy to be able to handle big hands where force may be applied by bariatric patients requiring the use of the handles for balance and stability
- Polyurethane wheels of easy moving.

Specifications

Weight Capacity	400 Kg
Width:	60 cm
Height of pan:	3.5-4 cm
Length:	60 cm
Lip Height:	.5 cm

Colour:

Warranty: 1 year limited