



# Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Raised Toilet Seating with Armrests



The raised seating Toilet Aid with armrests is designed to assist patients in lowering themselves onto the toilet, and raising from it.

### Features

- Features removable, rotating armrests and is height adjustable.
- Easy to clean.
- Fits most standard toilets.
- Secured at the rear of the seat using the existing holes on the toilet basin.
- Replaces existing seat and lid.
- Plastic moulded armrests.

### Specifications

<b>Weight Capacity</b>	114 kg
<b>Seat height:</b>	6cm, 10cm, 15cm
<b>Width:</b>	50 cm
<b>Length:</b>	475 cm
<b>Seat Width:</b>	36 cm
<b>Seat Depth</b>	42.5 cm
<b>Weight:</b>	4.3 Kg

**Warranty:** 1 year limited