



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: reception@mobileindependence.com.au

Web: www.mobileindependence.com.au

Reacher's



This lightweight Reacher has an easy to use trigger handle making it ideal for use around the home.

Features

- Magnetic head which allows user to pick up small metal objects such as pins
- Hook on the end of the Reacher which aids in dressing
- Palm lever which helps clamp the jaws for those who have poor hand strength
- Handy clip to attach to a rollator or wheelchair
- Rotating head to reach awkward items

Specifications

	Small	Medium	Large
Overall Length	60 cm	75cm	90cm
Weight	150gm	200gm	250gm

Warranty: 1 year limited