

Mobile Independence Pty Ltd

ABN: 90 531 747 431



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

**1/60 Kremzow Rd
Brendale Qld 4500**

Ph: 07 3881 3868

Fax: 07 3881 0525

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

SafeCare Floor Bed



The SafeCare Floor Bed has a custom design that allows the bed to lower to 75mm from the floor. Newly designed moulded columns at the ends of the bed are used for height adjustment and tilt.

Features

- Electric height and back adjustment
- Electric adjustable Trendelenburg tilt
- Personal hand control
- Head board colour Antique White
- Head and foot board (removable)
- Full clearance for lifting machines
- Maximum height 750 mm
- Minimum height 75 mm
- Mattress base 820 x 1920 mm

Options

- Electric adjustable knee bend
- Battery backup
- Patient lockout
- CPR button
- Self Help pole
- IV pole and holders
- Mountings for side rail use (free) must be ordered prior to manufacture
- Side rails are available but not necessary for this bed
- Wood grain head and foot boards

Warranty: 10 year Structural and 5 year Electrical