

Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Shower Chair



This Shower Chair is a height adjustable chair with a one-piece moulded plastic seat. It is designed to be placed in the shower (or bath) to provide safety and comfort while showering.

Features

- Height Adjustable
- Zinc Pre-treatment (anti-corrosive)
- One piece blow moulded seat
- Drainage holes in seat to help reduce the water build up
- Micro adjustments on legs for uneven floors
- Supplied fully assembled

Specifications

Weight Capacity	130 kg
Seat height:	41 - 53 cm
Width:	51 cm
Length:	48 cm
Seat Width:	44 cm
Weight:	5 Kg

Warranty: 1 year limited