



Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Side Folding Walker



This walker folds sideways, 8" Wheel.

Features

- Aluminium powder coated frame
- Compact & Easy to fold
- Multi height adjustable handles
- Large soft foldable seat storage bag
- Reversible & comfortable folding back rest
- Easy adjustable hand brakes
- Includes Cane Holder

Specifications

Weight Capacity	150 kg
Seat height:	58.4 cm
Folded Width:	23 cm
Weight:	6.4 kg

Colour: Red

Warranty: 1 year limited