

Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

Standard Hi Lo Bed



The standard Hi Lo is an aged care bed featuring cross Lift action to ensure the raising mechanism is within the horizontal plane, reducing the bed footprint and the possibility of damage to walls and other capital equipment.

Features

- Electric adjustable bed height
- Electric adjustable back rest
- Personal push button hand control
- High quality 'Dewert' German made electronics
- Removable head and foot boards
- Queen Anne head boards or choose from our range of shapes
- Self Help pole mounting bracket
- Wall bumper to keep bed distanced away from the wall whilst in use
- Hoist clearance 150 mm
- Mattress platform measurement: width: 900 mm x length: 1940 mm

Options

- Electric knee bend and foot raise
- Carer's control panel (foot end)
- Patient lock-out magnetic swipe key
- Battery back-up
- Under bed lighting
- Central locking castors
- Custom sized beds made to order
- Bed sticks - side attachable
- Pull out extension
- Mattress and / or mattress bolster in-fill for pull out extension
- Side rails (cot sides)
- Set of padded covers for side rails (cot sides)
- CPR button
- IV pole / Self Help pole
- Optional colours from our range for headboards
- Electric adjustable Trendelenburg tilt

Warranty: 10 year Structural and 5 year Electrical