



Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

T Shape Canes



Cyclone
Blue



Cyclone
Green



Rose



Blue
Swirl



Paisley

These T Shape Handle Canes are made of aluminium with a push button height adjustment and supplied with a soft grip, long lasting rubber handle. The handle provides a solid secure feel that conforms to the shape of your hand. These canes are available in Black, Bronze and Chrome.

Features

- Height adjustable in 25mm (1") increments
- Constructed with 22mm anodised aluminium
- Secured with locking ring for added safety
- Non slip rubber tips (with metal inserts) ensure a secure grip

Specifications

Weight Capacity 114 kg
Height: 74 - 97 cm
Weight: 370 g

Colour: Black, Bronze, Chrome, Cyclone blue, Cyclone Green, Rose, Blue Swirl, Paisley

Warranty: 1 year limited