

# Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Tri Walker 8" Wheels



4 wheeled, Aluminium, Loop Brake, 8" Wheel, Blue Patterned.

### Features

- Lightweight aluminium frame
- Four 8" x 1" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard

### Specifications

<b>Weight Capacity</b>	114 kg
<b>Width:</b>	68 cm
<b>Height:</b>	86 -99cm
<b>Length:</b>	58 cm
<b>Weight:</b>	4.9 kg

**Colour:** Red/Blue Patterned

**Warranty:** 1 year limited