



Mobile Independence Pty Ltd

ABN: 90 531 747 431

3/16 Commercial Drive
Caboolture, QLD 4510

Ph: 07 3881 3868



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2015

REGN. Number 3553

Email: reception@mobileindependence.com.au

Web: www.mobileindependence.com.au

Venus



The Heartway Venus is a comfortable and easy to use mobility scooter, with a sleek exterior that will make you the envy of all your neighbours. This mobility scooter is fitted with an easy to read and easy to operate control panel that displays the speedo, clock, climate indicator, odometer, tripmeter and battery charge indicator. You will always be up to date when you're driving your Venus mobility scooter. The Venus also features double wishbone, independent front and rear suspension system that provides a high level of comfort for the luxurious mobility scooter driver.

Specifications

Weight Capacity	135 kg
Width:	68 cm
Height:	132 cm
Length:	127 cm
Seat Width:	45 cm
Ground Clearance:	10 cm
Weight:	105 kg