

# Mobile Independence Pty Ltd

ABN: 90 531 747 431



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

**1/60 Kremzow Rd  
Brendale Qld 4500**

**Ph: 07 3881 3868  
Fax: 07 3881 0525**

**Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)**

**Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)**

## Ward Bed



The Ward bed is an aged care bed featuring cross Lift action to ensure the raising mechanism is within the horizontal plane, reducing the bed footprint and the possibility of damage to walls and other capital equipment.

### Features

- Electric adjustable bed height
- Electric adjustable back rest
- Electric adjustable Trendelenburg tilt
- Personal push button hand control
- Patient lock-out magnetic swipe key
- Removable head boards
- Electric knee bend and foot raise
- CPR emergency button
- Self Help and IV pole mounting bracket
- Set of side rails (cot sides)
- Central locking twin castors (125 mm)
- Battery backup
- High quality 'Dewert' electronics
- Round corner buffers for all four corners
- Hoist clearance 150 mm

### Options

- Padded covers for head boards and/or side rails (cot sides)
- Side attachable bed sticks
- IV Pole
- Self Help pole
- Head and foot boards are also available in various shapes with wood grain or coloured laminate finishes
- Extra Long 'fixed' length mattress platform OR a 250 mm pull out mattress platform extension
- Mattress and/or mattress bolster in-fill for pull out extension
- Under bed lighting
- 150 mm central locking castors to allow for Gazunder use
- Upgrade to SWL 350 kg King Single/Custom size• Electric adjustable Trendelenburg tilt

**Warranty: 10 year Structural and 5 year Electrical**