



# Mobile Independence Pty Ltd

ABN: 90 531 747 431

3/16 Commercial Drive  
Caboolture, QLD 4510

Ph: 07 3881 3868

Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## 20", 22" & 24" Wheelchair



Suitable for hospitals, nursing homes, rental use.

### Features

- Heavy Duty Construction
- Spill proof upholstery
- Solid tyres
- Flip up removable armrest
- Detachable swing away leg rest
- Wide front castors

### Specifications

<b>Weight Capacity</b>	250 kg
<b>Seat height:</b>	52 cm
<b>Width:</b>	75 cm
<b>Height:</b>	97-108 cm
<b>Length:</b>	106 cm
<b>Seat Width:</b>	50, 56, 61 cm
<b>Seat Depth:</b>	44 cm
<b>Weight:</b>	23.3 kg

**Colour:** Grey

**Warranty:** 1 year limited