



Mobile Independence Pty Ltd

ABN: 90 531 747 431



1/60 Kremzow Rd
Brendale Qld 4500

Ph: 07 3881 3868
Fax: 07 3881 0525

QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: sales@mobileindependence.com.au

Web: www.mobileindependence.com.au

20", 22" & 24" Wheelchair



Suitable for hospitals, nursing homes, rental use.

Features

- Heavy Duty Construction
- Spill proof upholstery
- Solid tyres
- Flip up removable armrest
- Detachable swing away leg rest
- Wide front castors

Specifications

Weight Capacity	250 kg
Seat height:	52 cm
Width:	75 cm
Height:	97-108 cm
Length:	106 cm
Seat Width:	50, 56, 61 cm
Seat Depth:	44 cm
Weight:	23.3 kg

Colour: Grey

Warranty: 1 year limited