



# Mobile Independence Pty Ltd

ABN: 90 531 747 431

1/60 Kremzow Rd  
Brendale Qld 4500

Ph: 07 3881 3868  
Fax: 07 3881 0525



QUALITY ASSURED COMPANY BY AS/NZS ISO 9001:2008

REGN. Number 3553

Email: [sales@mobileindependence.com.au](mailto:sales@mobileindependence.com.au)

Web: [www.mobileindependence.com.au](http://www.mobileindependence.com.au)

## Wombat Lifter



The Wombat is a compact hoist for home or institutional use. It is made to the same durable standard as the other lifters in our range. Electric lift and foot operated leg opening

### Features

- Tilt cradle frame. This provides maximum comfort and safety for patient and carer.
- Smaller front castors can be supplied when under bed clearance is an issue.
- Standard 100mm castors are supplied as standard for easy rolling.
- Haycomp hoists are designed to pick up from the floor and lift to clear 1 metre. The clearance height may vary with different hanger and sling combinations.
- Australian Made

### Specifications

<b>Weight Capacity</b>	140 Kg
<b>Lift:</b>	Electric
<b>Leg Opening:</b>	Foot Pedal
<b>Front Leg Height:</b>	14 cm
<b>Base Width:</b>	65 cm
<b>Overall Length:</b>	98 cm
<b>Mast Height:</b>	135 cm
<b>Floor Lift:</b>	Yes

### Specifications

Power Tilt

It can be supplied with a yoke hanger for use with GP slings

Weigh scale

Quick release version that can be easily broken down to enable transport in a car.

**Warranty:** 2 year limited